## LANTOS BY THE NUMBERS

**OUR VISION:** TO REPLACE EARMOLD IMPRESSIONS WITH SCANNING TECHNOLOGY



22,000 PEOPLE SCANNED

## BY THE LANTOS 3D EAR SCANNING SYSTEM

across a diverse set of age groups, geographies, and hearing conditions without a single adverse event.

> Lantos 3D Ear Scanning System is the only FDA 510K cleared and registered ear scanning device.

3,100,000 Earmold impressions taken globally

each year.

people prefer scanning to earmold impressions


**4**,700



from Lantos scans for hearing aids, professional audio, consumer audio and hearing protection.

## **1,000,000+** DATA POINTS

captured per ear by the Lantos 3D Ear Scanning System. Individual images are captured and converted through a proprietary algorithm to generate an accurate real-time 3D model of the full ear.



## ONE SCAN. ENDLESS DEVICES.

A single Lantos Scan can generate a truly limitless array of custom hearing devices.



2019

**DIGITAL EAR IMPRESSIONS** BREAKTHE MOLD.



